

Dewatering Device



- **Squeeze filtrate out of high-solids pulp**
- **Standardize sample preparation for pulp testing**
- **Stabilize hot stock samples by removing excess chemicals**
- **Convenient sample size “cakes” for storage, shipping and handling**

Description

The Pulmac Dewatering Device is a tool designed to allow the controlled dewatering of an unknown (3-15% concentration) pulp sample to a consistent solids content in just five minutes. Weighing out pulp from this dewatered sample accurately “predicts” an oven dried equivalent pulp sample for a variety of test purposes.

Pulmac International
A Division of Cowan Technologies Inc.
32 East State Street, Suite 10
Montpelier, Vermont
05602 USA

Phone: + 1.802.223.2663
Fax: + 1.802.223.2664
Email: info@pulmac.us <<mailto:info@pulmac.us>>
URL: www.pulmac.us <<http://www.pulmac.us>>

Represented by: REGMED Indústria Técnica de Precisão
Av. Dr. Alberto Jackson Byington, 1595
CEP 06276-000, Osaco - Brazil
Phone + 55 11 3601-5700
Fax. + 55 11 3601-5973
regmed@regmed.com.br
www.regmed.com.br