



A Division of Cowan Technologies Inc.

# Accept Collection Basket



Pulmac Accept Collection Basket

**The Pulmac Accept Collection Basket is used to collect fiber accepted from the Pulmac MasterScreen.**

## Description

ACB (Accept collection basket) captures all fiber that will not pass a 200 mesh screen. Separate frame and basket. Will capture up to 100 grams, oven dried equivalent, of pulp from the MasterScreen. Basket can be taken out of frame and backwashed to remove unwanted fiber or to flush fiber into another container.



Represented by:

## Benefits

Screened fiber can be captured for performing kappa tests, making handsheets, or examining debris in accepted fraction of pulp sample. can be calculated by subtracting the fraction it.

## Space Required:

- Dimensions: 60 x 32 x 42 cm ( W x L x H )
- Weight: 12 kg

Pulmac International  
A Division of Cowan Technologies Inc.  
32 East State Street, Suite 10  
Montpelier, Vermont  
05602 USA

Phone: + 1.802.223.2663  
Fax: + 1.802.223.2664  
Email: [info@pulmac.us](mailto:info@pulmac.us)  
URL: [www.pulmac.us](http://www.pulmac.us)

Represented by:

REGMED Indústria Técnica de Precisão Ltda.  
Av. Dr. Alberto Jackson Byington, 1595  
CEP 06276-000  
Osasco - Brazil  
Phone + 55 11 3601-5700  
Fax + 55 11 3601-3973  
[regmed@regmed.com.br](mailto:regmed@regmed.com.br)  
[www.regmed.com.br](http://www.regmed.com.br)